DENSORAPID D

DENSORAPID D

Transmission Densitometer according to ISO 5, ISO 14807, ASTM E1079



For precise determination of optical densities of X-ray films

DENSORAPID D Transmission Densitometer, hand held

For precise determination of optical densities of X-ray films the DENSORAPID D is equipped with innovative digital measuring technique.

The illuminated display enables easy readings even in darkened rooms and calibration of the instrument, done with the software, is simple and much more accurate than on conventional densitometers.

Fault-free measurements even in slightly inclined positions are guaranteed by the ergonomic hand probe with a soft tip avoiding scratches on the film.

The battery is charged with a standard USB charger and ensures an operating time of more than 100 hours with display illumination switched off.

Specification:

Housing:	ABS plastics, IP 65
Hand probe:	Aluminum, anodized
Measuring range:	D 0.00 > 5.00
Resolution:	D 0.01
Accuracy:	+/- D 0.03
Repeatability:	+/- D 0.02
Aperture size:	Ø 2.0 mm
Communication port:	USB 2.0, Micro USB
Battery operating time:	100 hours (illumination off)
USB Charger:	100 240 V, 50 60 Hz, Micro USB
Dimensions unit:	135 x 73 x 29 mm (LxWxH)
Dimensions probe:	Ø 14 x 120 mm (ØxL), cable 1 m
Dimensions case: :	325 x 270 x 55 mm
Weight unit:	0.25 kg
Weight total with case:	0.90 kg
weight total with case.	0.90 KY

Required measuring light source - Film viewer:Color temperature:ca. 2,800 ... 5,000 KLuminance:ca. 70,000 ... 100,000 cd/m²

Scope of supply:

DENSORAPID D Densitometer complete with: USB Charger Transport case Measuring mask, magnetic, 10x24 cm / Ø 10 mm KOWOSTEP D Density step tablet 7 density values, >0,30 ... >5,00 Re-calibration interval 3 years Test certificates, Certificates of Conformity

Accessories:

Measuring mask, magnetic, 10x24 cm / Ø 10 mm

X Densitometer holder for KOWOLUX X universal for all type hand held Densitometers

BAM X-001 Density step tablet 15 density values, >0,20 ... >5,00 Re-calibration interval 5 years









PKT_1340540-2203-03_EN

Article No.: 13 40084 Article No.: 13 40087

Article No.: 13 40540

Article No.: 13 40546

Article No.: 13 40547

Article No.: 13 40104 Article No.: 13 40207

Article No.: 13 40205



Solinger Strasse 186 40764 Langenfeld +49 2173 22383 Römerstrasse 12 83308 Trostberg +49 8621 979110 Germany www.kowotest.de info@kowotest.de